Deta
Elis
Holding
Gesundheit
Innovationen
Mensch

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A Skin Care Device Based On
Quantum Effect Combined With
Harmonization And Correction Of
Natural Electromagnetic Emission
Of Natural Biosystems

DeVita
Cosmo-50

Introduction

Purpose of the Device: The DeVita Cosmo-50 is a generator of weak electromagnetic radiation, identical to that of natural organic systems, as well as optical radiation. It is intended for skin care based on the harmonization and correction of natural electromagnetic radiation of natural biosystems.

Principle of Action: The device's corrective action utilizes frequencies based on the methods of R. Rife, R. Voll, and F. Morrel. It is designed for individual use, in cosmetic clinics and at home.

Important Safety Information

Warning: The device may only be used by children under adult supervision. Children are not allowed to play with the device. Persons with physical or mental disabilities may use the device only under the supervision of appropriate personnel.

Safety Profile: The device is non-toxic, non-flammable, non-explosive, and does not emit harmful radiation.

Package Contents

The package contains the following items:

"DeVita COSMO" device.....	1 pc.
Remote QUANT unit.....	1 pc.
Power elements (AA batteries, 1.5V).....	2 pcs.
Charging cable for the QUANT unit: USB A to Type C.....	1 pc.
USB Type-C to Type-C cable.....	1 pc.
Passport/Operating Manual	1 pc.
Guidelines.....	1 pc.
Packaging.....	1 pc.

Technical Specifications

1. Maximum number of programs, units.....	50
2. Frequency range.....	0,1÷1000Hz
3. Optical range (QUANT Unit).....	468, 525, 590, 660, 880nm
4. Continuous Operation Time, not less than.....	30 hours
5. Display type.....	TFT 128x160
6. Power (Main Device).....	2 x AA Batteries(1.5V)
7. Main Device Dimensions.....	135 x 70 x 25 mm - Weight: ≤ 0.2 kg

8. QUANT Unit Dimensions.....95 x 47 x 27 mm - Weight: ≤ 0.05 kg

Environmental protection

Disposal: The built-in battery contains environmentally hazardous substances. At the end of the device's life, the battery must be removed by a service center before recycling. 

Maintenance, Storage, and Transportation

The device without packaging should be stored indoors at 10°C to 35°C with relative humidity $\leq 80\%$. Protect from moisture, aggressive chemicals, and mechanical impact (drops/shocks). Clean the exterior with a soft, dry cloth. It is recommended to keep the packing container during the warranty period.

Warning! Do not disassemble the device.

Overview

Main Device Layout

1. Power "ON/OFF" button
2. Display
3. External 5V Power Input Socket (Note: For external power supply, not for charging AA batteries)
4. Navigation buttons (Up/Down)
5. Option selection button (OK)
6. The "RETURN" button
7. Operation indicator
8. Remote QUANT unit


Fig. 1

Operating

To Turn On, press the "ON/OFF" button (Fig. 1).


The display stays on for 30 seconds after the last key press. To turn the screen off manually, press the "ON/OFF" button. Press it again to wake the screen.

Auto-Shutoff: To save power, the screen turns off automatically after 10 seconds of inactivity (configurable).

After pressing the "ON/OFF" button (Fig. 1), a splash screen will appear on the display. In 2 seconds, the device will automatically enter the operational mode set during the last shutdown.

Navigating the main "Menu"


Submenu «Settings»

  Use the navigation buttons to select "Settings".

 Press the "OK" button; The "Settings" submenu will appear on the display.

Language selection

  Use the navigation buttons to select "Language".

 Press the "OK" button.

The "Language" submenu will appear on the display.


  Use the navigation buttons to choose the desired language.

 Press the "OK" button.

 Press "RETURN" to navigate back to the "Settings" submenu.


Sound setting

  Use the navigation buttons to select "Sound".

 Press the "OK" button.

The "Sound" submenu will appear on the display.


  Use the navigation buttons to select "On/Off".

 Press the "OK" button.

 Press the "RETURN" button to navigate back to the Settings submenu.


Note: Selecting "Off" disables also the keyboard sounds.

About the device

  Use the navigation buttons to select "About the Device".

 Press the "OK" button.

The logo will appear on the display, followed by the model name, serial number and the software version of your DeVita Cosmo device.

QUANT connection

Ensure that the QUANT unit number matches the device number.

This connection procedure is necessary if the Quantum Head has been disconnected, damaged or lost and needs to be replaced with a new one.

Before connecting the Quantum Head, you must first charge it for at least two hours, using the USB Type-A to Type-C cable.

Make sure that the batteries in the DeVita COSMO device must not be discharged.

70% 
Menu
Programs
Settings

70% 
Settings
Language
Sound
About the device
QUANT connection


To enable the connection of the Quantum Head to the COSMO device, use the USB Type-C to Type-C cable.

To program, follow these steps:

1. Turn on the device
2. Enter the "Menu" window by pressing the Back key
3. Select "Settings" and press OK
4. Select "Connecting QUANT" and press OK
5. Connect the Quantum Head to the COSMO device using the USB Type-C to Type-C cable.
6. Note that the LEDs on the Quantum Head should light up – indicating programming mode.

- If the LEDs do not light up, rotate the USB connector on the device 180 degrees.
- If the LEDs still do not light up, reverse the USB connector on the Quantum Head as well.

Remember or mark the positions of the USB connectors for future programming.


Once you have entered programming mode:

7. To connect the head, select "Connect?" and press OK.

8. If you need to disconnect the head, select "Disconnect?" and press OK.


9. Note that the LEDs will briefly flash to confirm the programming selection.

Submenu "Device"

The device operates in two modes:

1. COSMO: Standalone operation.
2. COSMO+QUANT: Combined operation with the remote QUANT unit via Bluetooth.


Use the navigation buttons to select either "COSMO" or "COSMO+QUANT".

Press "OK" to proceed to the programs menu.

Note: In the combined mode, program execution time is automatically reduced by half due to the enhanced effect of combined electromagnetic and quantum effects

Program Selection


Use the navigation buttons to select the desired program.

Press the "OK" button.

The program then will start and the operation indicator (7 - fig.1) will flash green.

Follow the specific program's guidelines about placement and QUANT unit.

A submenu for the ongoing program will appear on the display.

A countdown timer of the program appears at the top of the display.


In order to lower battery drain, the device screen automatically turns off after 10 seconds of inactivity.

To connect, please ensure the QUANT unit is within 2 meters of the main device. When the Bluetooth connection is established with the QUANT unit, the display shows the battery icon, the charge level of the remote unit's battery and the remaining charge percentage, as well as the Bluetooth standard icon.

Troubleshooting: If the connection fails, the display will read "No connection with QUANT". Press "OK" to exit and try again.


Information on program execution

To obtain information on the program's progress while the screen is off, press the "ON/OFF" button (see Fig. 1) .

A submenu with the current program's status will appear on the display.


Stop / Pause the program

To stop the program,

  Use the navigation buttons to select "Yes".


 Press the "OK" button.

To pause the program during execution:

  Use the navigation buttons to select "Pause".

 Press the "OK" button.

A submenu with "Pause" will appear on the display. To quit from the pause:

  Use the navigation buttons to select "Yes".

 Press the "OK" button.


Note: The device automatically turns off if left paused for 30 minutes.

Device Placement

Immediately after starting the program, place the device with the back side against your body, (e.g., in your pocket). Then follow the specific program guidelines.

Power & Battery Management

The main device (DeVita COSMO) requires 2 x AA Batteries (1.5V) to operate. The display shows a battery icon with the charge level (white). On low battery the indicator turns red.


If completely discharged, the screen displays "Change the battery!" and the device turns off.

Replacing the batteries

Locate the battery compartment on the back of the device. Press the latch to open the cover, remove the depleted and insert new batteries.

Note: Stored programs are saved during battery replacement.

Battery discharge indicator for the remote QUANT unit

The battery charge indicator is lit in blue and appears as a battery icon showing the charge level and the remaining charge percentage.

If the battery charge is low, the indicator lights up red.

Battery charging

The remote QUANT unit has a built-in rechargeable battery. To charge, connect the USB Type-C end of the provided cable to the QUANT unit's charging port. Connect the USB-A end to a computer or a 5V/1A power adapter. The indicator lights blue when charged. It turns red when the battery is low.

Manufacturer's Warranty

The manufacturer guarantees that the "DeVita COSMO" device meets the technical specifications provided that the user complies with the rules of operation, transportation and storage. The warranty period of the device is 24 months from the date of purchase from a retail store.

If no proof of purchase is available, the period is calculated from the production date by the manufacturer.

During the warranty period, the owner is entitled to free repairs upon presentation of the warranty repair form. Warranty repairs are conducted at the manufacturer's facilities. Transportation of the defective device is at the expense of the buyer.

Claims are not accepted, and warranty repairs are not performed in case of damage to the device's body and/or without presenting an application for warranty repair.

Warranty does not cover damage caused by:

- defects due to improper use;
- defects caused by natural disasters;
- If the device seals are broken;
- If there is mechanical damage (cracks, chips, etc.).

The buyer is entitled to have the defective device replaced with a new one during the warranty period in the following cases:

- the device has been repaired three times during the warranty period;
- the device is unrepairable.

Warning! The manufacturer reserves the right to amend changes to the design of the device and software without prior notice, without compromising existing characteristics.

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	TECHNICAL FILE EC DECLARATION OF CONFORMITY
---	--

EC Declaration of Conformity Nr.W15

The company **DETA-ELIS EUROPA GMBH** with a seat at Kaiserstr. 74, 63065 Offenbach, Germany

Herewith to DECLARE

That the equipment described below complies with all applicable at the time of construction essential requirements of the **Directive 2011/65/EU (RoHS)** and the **Directive 2014/30/EU (EMC)**.

This declaration relates exclusively to the equipment in the state in which it was placed on the market, and excludes alterations or manipulations, components which are added and/or operations carried out subsequently by the final user without our official consent.

Equipment: DeVita
Type / Model: Cosmo
Standards Applied: EN 55014-1:2006+A1:2009+A2:2011
EN 55014-2:1997+A1:2001+A2:2008
EN 61000-3-2:2006+A1:2009+A2:2009
EN 61000-3-3:2008

Date: 23 January 2018

Authorized Representative: The General Manager Ms. Konopleva Tatyana

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CONCLUSION: IN REFER TO SGS TEST REPORT NO.: H0FA0004
NO DEVIATION FROM THE STANDARDS HAVE
BEEN IDENTIFIED. THEREFORE THE TESTED
PRODUCT IS IN COMPLIANCE WITH THE
APPROPRIATED STANDARDS.

This Verification of Compliance has been granted to the applicant based on the results of tests, performed by Laboratory of SGS-CSTC Standards Technical Services Co., Ltd. on sample of the above-mentioned product in accordance with the provisions of the relevant specific standards. The CE marking as shown below can be affixed, under the responsibility of the manufacturer, after completion of an EC Declaration of Conformity and compliance with all relevant EC Directives.

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company, signature 1 and 2

B. Ladd

SGS Germany GmbH, Hofmannstr. 50, D-81379 Munich is testing facility for
ELECTRO-MECH. COMPONENTS TESTING ENVIRONMENTAL SIMULATION ELECTROMAGNETIC COMPATIBILITY
MEDICAL DEVICE TESTING BATTERY TESTING PRODUCT SAFETY TELECOM CONFORMANCE TESTS

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Member of the SGS Group (Société Générale de Surveillance)

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Deta-Elis Europa GmbH